

Child & Adult
Nasal Speculum




info@vtspeculum.com




Rue du Bailli 59
1050 Ixelles, Belgium



www.vtspeculum.com



VISION TECH



Introducing LED Nasal Speculum from Vision Tech

Brand new speculum sets a new benchmark in ergonomics. Following the User-Centered Design, over quite a long period, LED Nasal Speculum gained a full understanding of the needs, wants and frustrations of end users. Subsequent painstaking attention to detail at every stage of the design process resulted in the creation of the ideal handles and the embedded light.

Precision and expert manufacture provides quality, accuracy and longevity. The addition of highly-praised design brings beauty and elegance.



Ergonomic Design



Light Weight



Anti-Slip Grips




Vision Assistance LED



1hr Usage Time



Rechargeable Battery



1hr Usage Time

Highly optimized circuitry and a high quality battery ensures a full 1 hour of usage time.

Light Weight

The perfect blend of stainless steel and rubber plastic that gives the user superior comfort also ensures the light weight.

Cleaning

There is no need to disassemble the VT Nasal Speculum after each use. Clean the assembled speculum like a regular speculum using a clinical agent.

Anti-Slip Grips

Anti-slip rubber surface provides a better grip. Small holes on the hard plastic cover provide an anti-slip grip and superior operator comfort.

Rechargeable Battery

A quick pull tab releases a charging port to recharge the battery. A Micro-USB port gives a universal port for charging and a small LED indicator light shows when the battery is fully charged.

Vision Assistance LED

The high intensity LED inside gives vision assistance and full control over where to point the light. As compared to wearing on the head, this light module brings out a clear cool white exposure at 5000 Kelvin.