

SimplyGo Mini

Portable Oxygen Concentrator

Sleep and Respiratory Care
SimplyGo Mini Launch Pack
June 2015

PHILIPS

SimplyGo Mini

At a glance

INTUITIVE

Easy to use and to read touch screen

OXYGEN DELIVERY

5 Pulse dose settings
1000 ml/min oxygen capacity

LIGHTWEIGHT DESIGN

Only 2,3kg (<5lbs) with standard battery

VERSATILITY

Two battery configurations:
- Standard
- Extended

RUGGED

Impact resistant design





SimplyGo Mini

Gives patients the
confidence to go
further

PHILIPS



SimplyGo Mini

Lightweight design

Today's oxygen patients are living their life, taking trips to the store or spending weekend away with friends and family.

SimplyGo Mini is a powerful choice for a lightweight portable oxygen concentrator for those patients.

With a weight of only 2,3kg* (5lbs) they now have an option to do what they want, when they want.

* With the standard battery

PHILIPS

SimplyGo Mini


Intuitive design

Easy to use touch screen interface to:

- Check battery status,
- Set desired settings
- Set language or screen brightness

No more unexplained/confusing alarms

In case of an alarm, SimplyGo Mini will show a written message to give patients explanations about the issue and offer troubleshooting options



SimplyGo Mini

Versatility

2 battery options available



	SIMPLYGO MINI		SIMPLYGO MINI
Weight (device + battery)	2,3kg (5,0 lbs)	Weight (device + battery)	2,7kg (6,0lbs)
Battery Autonomy (Setting 2)	Approx. 4,5 hrs	Battery Autonomy (Setting 2)	Approx. 9,0 hrs
Standard battery weight	0,5kg (1,0 lbs)	Standard battery weight	0,95 kg (2,0 lbs)

SimplyGo Mini accessories

Versatility



External battery charger



Standard battery



Extended battery



**Backpack
(Black or Brown)**



**Carry case
(Black or Brown)**



**Accessory Bag
(Black or Brown)**



SimplyGo Mini

Designed to be as
reliable as the
SimplyGo

PHILIPS

SimplyGo Mini

How can reliability benefit your business?



Increase in employee
value-added activities because
of fewer failure calls



Decrease in resources
and overtime payments
to deal with failed equipment
Calls and pick-ups



Decrease in inventory
Because of fewer
swap-outs



Decrease in resources
and overtime payments
to deal with failed equipment
Calls and pick-ups



SimplyGo Mini

Warranty

STANDARD WARRANTY

- 3 years on SimplyGo Mini
- 1 year on sieve-beds
- 90 days on battery

EXTENDED WARRANTIES

- 3 years on SimplyGo Mini – 3 years on the sieve-Beds
- 2 year sieve warranty extension

A white Philips SimplyGo Mini CPAP machine is shown in the foreground, resting on a dark surface. The machine has a clear tube connected to the top. The Philips logo and 'RESPIRONICS' are visible on the front. The background is blurred, showing a man and a woman in a domestic setting.

PHILIPS
RESPIRONICS

SimplyGo Mini

Two Packaging options



SIMPLYGO MINI – STANDARD BATTERY PACK

Pack will include: Includes SimplyGo Mini portable oxygen concentrator, carrying case, one **standard** rechargeable Lithium Ion battery, DC and AC power cords, accessory bag, and user manual



SIMPLYGO MINI – EXTENDED BATTERY PACK

Pack will include: Includes SimplyGo Mini portable oxygen concentrator, carrying case, one **extended** rechargeable Lithium Ion battery, DC and AC power cords, accessory bag, and user manual

SimplyGo & SimplyGo Mini

SimplyGo

Our smallest Continuous Flow Portable Concentrator

Oxygen delivery

- Pulse mode: 6 settings
- Continuous flow : 0,5 to 2LPM

Enables:

Day to day ambulatory usage:

- Less than 5 kg
- Easily transportable

Night usage, for weekends, vacation...

- Continuous flow
- Sleep mode
- FAA approved for airline travel

Proven reliability on the field

SimplyGo Mini

Our lightweight POC with the highest oxygen output

Oxygen delivery

- Pulse mode: 5 settings
- 1000ml/min oxygen capacity

Day to day ambulatory usage:

- Less than 2,5 kg – Easily portable
- Up to 9hrs battery autonomy*
- Optional backpack for hands-free use

Designed for reliability

** With extended battery, setting 2, 2BPM*

